

SFT2D3 Antibody



PACO12153

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SFT2D3

Uniprot

Q58719

Synonyms:

SFT2 domain containing 3; SFT2D3; MGC5391 ;

Immunogen:

Human SFT2D3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A