

SH3BGRL3 Antibody



PACO12180

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SH3BGRL3

Uniprot

Q9H299

Synonyms:

SH3 domain binding glutamic acid, rich protein like 3; SH3BGRL3; SH3BP-1; TIP-B1 ;

Immunogen:

Human SH3BGRL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A