

SH3RF1 Antibody



PACO12190

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SH3RF1

Uniprot

Q7Z6J0

Synonyms:

SH3 domain containing ring finger 1; SH3RF1; FLJ21602; KIAA1494; POSH; RNF142; SH3MD2 ;

Immunogen:

Human SH3RF1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A