

SHMT1 Antibody



PACO12200

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SHMT1

Uniprot

P34896

Synonyms:

serine hydroxymethyltransferase 1 (soluble); SHMT1; CSHMT; MGC15229; MGC24556; SHMT ;

Immunogen:

Human SHMT1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A