

SIPA1L2 Antibody



PACO12216

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SIPA1L2

Uniprot

Q9P2F8

Synonyms:

signal-induced proliferation-associated 1 like 2; SIPA1L2; FLJ23126; FLJ23632; KIAA1389; SPAL2 ;

Immunogen:

Human SIPA1L2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A