

SIRPB1 Antibody



PACO12218

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SIRPB1

Uniprot

Q5TFQ8

Synonyms:

signal-regulatory protein beta 1; SIRPB1; CD172b; DKFZp686A05192; FLJ26614; SIRP-BETA-1 ;

Immunogen:

Human SIRPB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A