

SLC10A3 Antibody



PACO12236

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SLC10A3

Uniprot

P09131

Synonyms:

solute carrier family 10 (sodium/bile acid, cotransporter family), member 3; SLC10A3; DXS253E; P3 ;

Immunogen:

Human SLC10A3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A