

SLC16A12 Antibody



PACO12241

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC16A12

Uniprot

Q6ZSM3

Synonyms:

solute carrier family 16, member 12 (monocarboxylic acid, transporter 12); SLC16A12; CJMG; DKFZp686E188; MCT12 ;

Immunogen:

Human SLC16A12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A