

# SLC18A1 Antibody



PACO12245

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC18A1

**Uniprot**

P54219

**Synonyms:**

solute carrier family 18 (vesicular monoamine), member 1; SLC18A1; CGAT; VAT1; VMAT1 ;

**Immunogen:**

Human SLC18A1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A