

# SLC19A3 Antibody



PACO12246

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC19A3

**Uniprot**

Q9BZV2

**Synonyms:**

solute carrier family 19, member 3; SLC19A3; THTR2 ;

**Immunogen:**

Human SLC19A3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A