

SLC22A12 Antibody



PACO12252

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC22A12

Uniprot

Q96S37

Synonyms:

solute carrier family 22 (organic anion/urate transporter), member 12; SLC22A12; OAT4L; RST; URAT1 ;

Immunogen:

Human SLC22A12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A