

SLC22A2 Antibody



PACO12255

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC22A2

Uniprot

O15244

Synonyms:

solute carrier family 22 (organic cation transporter), member 2; SLC22A2; MGC32628; OCT2 ;

Immunogen:

Human SLC22A2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A