

SLC22A5 Antibody



PACO12256

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC22A5

Uniprot

O76082

Synonyms:

solute carrier family 22 (organic cation/carnitine transporter), member 5; SLC22A5; CDSP; FLJ46769; OCTN2; OCTN2VT ;

Immunogen:

Human SLC22A5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A