

SLC8B1 Antibody



PACO12258

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC8B1

Uniprot

Q6J4K2

Synonyms:

solute carrier family 24 (sodium/potassium/calcium exchanger), member 6; SLC24A6; FLJ22233; NCKX6; NCLX ;

Immunogen:

Human SLC24A6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A