

SLC25A10 Antibody



PACO12260

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC25A10

Uniprot

Q9UBX3

Synonyms:

solute carrier family 25 (mitochondrial carrier; dicarboxylate transporter), member 10; SLC25A10; DIC ;

Immunogen:

Human SLC25A10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A