

SLC25A20 Antibody



PACO12265

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

SLC25A20

Uniprot

O43772

Synonyms:

solute carrier family 25 (carnitine/acylcarnitine translocase), member 20; SLC25A20; CAC; CACT ;

Immunogen:

Human SLC25A20.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A