

SLC25A34 Antibody



PACO12270

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC25A34

Uniprot

Q6PIV7

Synonyms:

solute carrier family 25, member 34; SLC25A34; DKFZp781A10161; RP11-169K16.2 ;

Immunogen:

Human SLC25A34.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A