

# SLC25A6 Antibody



PACO12274

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

SLC25A6

**Uniprot**

P12236

**Synonyms:**

solute carrier family 25 (mitochondrial carrier; adenine nucleotide translocator), member 6; SLC25A6; AAC3; ANT3; ANT3Y; MGC17525 ;

**Immunogen:**

Human SLC25A6.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A