

SLC27A2 Antibody



PACO12276

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SLC27A2

Uniprot

O14975

Synonyms:

solute carrier family 27 (fatty acid, transporter), member 2; SLC27A2; ACSVL1; FACVL1; FATP2; HsT17226; VLACS; VLCS; hFACVL1 ;

Immunogen:

Human SLC27A2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A