

# SLC27A3 Antibody



PACO12277

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

SLC27A3

**Uniprot**

Q5K4L6

**Synonyms:**

solute carrier family 27 (fatty acid, transporter), member 3; SLC27A3; ACSVL3; FATP3; MGC4365; VLCS-3 ;

**Immunogen:**

Human SLC27A3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A