

SLC30A3 Antibody



PACO12280

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SLC30A3

Uniprot

Q99726

Synonyms:

solute carrier family 30 (zinc transporter), member 3; SLC30A3; ZNT3 ;

Immunogen:

Human SLC30A3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A