

SLC35A1 Antibody



PACO12284

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SLC35A1

Uniprot

P78382

Synonyms:

solute carrier family 35 (CMP-sialic acid, transporter), member A1; SLC35A1; CMPST; CST; hCST; inactive ;

Immunogen:

Human SLC35A1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A