

SLC37A4 Antibody



PACO12286

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SLC37A4

Uniprot

O43826

Synonyms:

solute carrier family 37 (glucose-6-phosphate transporter), member 4; SLC37A4; G6PT1; G6PT2; G6PT3; GSD1b; GSD1c; GSD1d; MGC15729; PRO0685; TRG19 ; ; G6PT

Immunogen:

Human SLC37A4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A