

SLC38A1 Antibody



PACO12287

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

SLC38A1

Uniprot

Q9H2H9

Synonyms:

solute carrier family 38, member 1; SLC38A1; ATA1; NAT2; SAT1; SNAT1 ;

Immunogen:

Human SLC38A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A