

# SLC39A4 Antibody



PACO12290

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC39A4

**Uniprot**

Q6P5W5

**Synonyms:**

solute carrier family 39 (zinc transporter), member 4; SLC39A4; AEZ; FLJ20327; MGC74741; ZIP4 ;

**Immunogen:**

Human SLC39A4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A