

SLC39A6 Antibody



PACO12291

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SLC39A6

Uniprot

Q13433

Synonyms:

solute carrier family 39 (zinc transporter), member 6; SLC39A6; LIV-1 ;

Immunogen:

Human SLC39A6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A