

SLC39A7 Antibody



PACO12292

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

SLC39A7

Uniprot

Q92504

Synonyms:

solute carrier family 39 (zinc transporter), member 7; SLC39A7; D6S115E; D6S2244E; H2-KE4; HKE4; KE4; RING5; ZIP7 ;

Immunogen:

Human SLC39A7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A