

SLC44A1 Antibody



PACO12297

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC44A1

Uniprot

Q8WWI5

Synonyms:

solute carrier family 44, member 1; SLC44A1; CD92; CDW92; CHTL1; CTL1; RP11-287A8.1 ;

Immunogen:

Human SLC44A1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A