

SLFN1 Antibody



PACO12313

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

SLFN1

Uniprot

Q499Z3

Synonyms:

schlafen-like 1; SLFN1; FLJ23878; MGC43873 ;

Immunogen:

Human SLFN1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A