

SLITRK3 Antibody



PACO12314

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLITRK3

Uniprot

O94933

Synonyms:

SLIT and NTRK-like family, member 3; SLITRK3; KIAA0848; MGC138681 ;

Immunogen:

Human SLITRK3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A