

PACO12318

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLN

**Uniprot**

O00631

**Synonyms:**

sarcosin; SLN; MGC12301; MGC125854; MGC125855 ;

**Immunogen:**

Human SLN.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A