

# SNCA Antibody



PACO12361

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SNCA

**Uniprot**

P37840

**Synonyms:**

synuclein, alpha (non A4 component of amyloid precursor); SNCA; MGC110988; NACP; PARK1; PARK4; PD1 ; Alpha-synuclein

**Immunogen:**

Human SNCA.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A