

SND1 Antibody



PACO12363

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SND1

Uniprot

Q7KZF4

Synonyms:

staphylococcal nuclease and tudor domain containing 1; SND1; TDRD11; p100 ;

Immunogen:

Human SND1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A