

PACO12374

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SNRPN

Uniprot

P63162

Synonyms:

small nuclear ribonucleoprotein polypeptide N; SNRPN; DKFZp686C0927; DKFZp686M12165; DKFZp761I1912; DKFZp762N022; FLJ33569; FLJ36996; FLJ39265; HCERN3; MGC29886; PWCR; RT-LI; SM-D; SMN; SNRNP-N; SNURF-SNRPN ;

Immunogen:

Human SNRPN.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A