

SNX15 Antibody



PACO12377

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SNX15

Uniprot

Q9NRS6

Synonyms:

sorting nexin 15; SNX15; HSAF001435 ;

Immunogen:

Human SNX15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A