

SNX20 Antibody



PACO12380

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SNX20

Uniprot

Q7Z614

Synonyms:

sorting nexin 20; SNX20; MGC35578; SLIC-1; SLIC1 ;

Immunogen:

Human SNX20.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A