

SNX27 Antibody



PACO12382

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SNX27

Uniprot

Q96L92

Synonyms:

sorting nexin family member 27; SNX27; KIAA0488; MGC126871; MGC126873; MGC20471; MRT1; MY014 ;

Immunogen:

Human SNX27.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A