

SORCS2 Antibody



PACO12394

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SORCS2

Uniprot

Q96PQ0

Synonyms:

sortilin-related VPS10 domain containing receptor 2;

Immunogen:

Human SORCS2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A