

# SOX13 Antibody



PACO12399

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SOX13

**Uniprot**

Q9UN79

**Synonyms:**

SRY (sex determining region Y)-box 13; SOX13; ICA12; MGC117216; Sox-13 ;

**Immunogen:**

Human SOX13.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A