

# SP6 Antibody



PACO12407

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SP6

**Uniprot**

Q3SY56

**Synonyms:**

Sp6 transcription factor; SP6; EPFN; EPIPROFIN; KLF14; MGC119662; MGC119663; MGC119664; MGC119665 ;

**Immunogen:**

Human SP6.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A