

SPA17 Antibody



PACO12408

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SPA17

Uniprot

Q15506

Synonyms:

sperm autoantigenic protein 17; SPA17; CT22; SP17; SP17-1 ;

Immunogen:

Human SPA17.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A