

SPAG8 Antibody



PACO12416

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SPAG8

Uniprot

Q99932

Synonyms:

sperm associated antigen 8; SPAG8; BS-84; HSD-1; MGC26201; SMP1; SPAG3; hSMP-1 ;

Immunogen:

Human SPAG8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A