

SPANXC Antibody



PACO12418

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SPANXC

Uniprot

Q9NY87

Synonyms:

SPANX family, member C; SPANXC; C; CT11.3; CTp11; SPANX-C ;

Immunogen:

Human SPANXC.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A