

# SPATA19 Antibody



PACO12423

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SPATA19

**Uniprot**

Q7Z5L4

**Synonyms:**

spermatogenesis associated 19; SPATA19; FLJ25851; SPAS1; spergen1 ;

**Immunogen:**

Human SPATA19.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A