

SPTBN2 Antibody



PACO12445

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SPTBN2

Uniprot

O15020

Synonyms:

spectrin, beta, non-erythrocytic 2; SPTBN2; GTRAP41; SCA5 ;

Immunogen:

Human SPTBN2.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A