

SPTLC2 Antibody



PACO12447

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SPTLC2

Uniprot

O15270

Synonyms:

serine palmitoyltransferase, long chain base subunit 2; SPTLC2; KIAA0526; LCB2; SPT2 ;

Immunogen:

Human SPTLC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A