

PACO12455

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SRF

Uniprot

P11831

Synonyms:

serum response factor (c-fos serum response element-binding transcription factor); SRF; MCM1 ;

Immunogen:

Human SRF.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A