

SRFBP1 Antibody



PACO12456

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SRFBP1

Uniprot

Q8NEF9

Synonyms:

serum response factor binding protein 1; SRFBP1; BUD22; FLJ25286; P49; Rib1; STRAP; p49/STRAP ;

Immunogen:

Human SRFBP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A