

SRP14 Antibody



PACO12461

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SRP14

Uniprot

P37108

Synonyms:

signal recognition particle 14kDa (homologous Alu RNA binding protein); SRP14; ALURBP; MGC14326 ;

Immunogen:

Human SRP14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A