

# SS18L2 Antibody



PACO12471

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SS18L2

**Uniprot**

Q9UHA2

**Synonyms:**

synovial sarcoma translocation gene on chromosome 18-like 2; SS18L2; KIAA-iso ;

**Immunogen:**

Human SS18L2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A